



Curva Lounger



Shown with optional sliding shelf and SoVi System

Perfect for:

- VibroAcoustic Therapies
- IV Drip Therapies
- Waiting Areas
- Salt Rooms


The **Curva Lounger** creates the perfect balance between comfort and sociability. The NASA inspired curve and tilt creates an incredibly comfortable position that is enhanced by our cloud soft 4" AeroCel padding system. But knowing that many times people want to interact with their friends we made sure that people are upright enough to make it easy to converse with other people.



SoVi Vibroacoustic System
(Optional)



Color LED Lighting
(Optional)



Integrated Heated Top
(Optional)

SPECIFICATIONS

LOUNGER SPECIFICATIONS

Top	Replaceable 4" Top Pad with Neck Bolster
Seat Height	18" (46 cm)
Weight Capacity	500 lbs. (227 kg)
Fabric	TerraTouch - Medical Grade, ISO 10993
Finish	Laminate
Padding	4" (10 cm) AeroCel™
Width	31" (79 cm)
Length	Pad length 75" (109 cm)
Table Weight	125 lbs.
Options (must be ordered with table)	SoViVibroacoustic System, Heated Top, Sliding Shelf, Power Outlet, Colored LED Lighting
Warranty	3 years parts, 2 years labor

LAMINATE COLORS




TERRATOUGH FABRIC COLORS



REPLACEABLE PAD


The pad is easily replaced if it ever gets damaged.



SOVI VIBROACOUSTIC SYSTEM

(OPTIONAL)


The system is designed to deliver all the latest forms of healing modalities via a very flexible integrated delivery system.



HEATED TOP

(OPTIONAL)


The same system we use on our top of the line spa tables, our system radiates an even heat over the entire body.



SLIDING SHELF

(OPTIONAL)


Designed to be stored out of the way or pulled out to either side, this acrylic shelf is easy to clean.



COLORED LED LIGHTING

(OPTIONAL)

A wide range of colors and patterns are available and are accentuated by the optional shelf.



EXTERNAL OUTLET

(OPTIONAL)

External outlet with USB port for charging devices.



Contact Oakworks
717.235.6807
www.shop.spatables.com
customerservice@oakworks.com